

Lens Tissues

- Soft, Premium-Grade Sheets
- 5 (MC-5) or 50 (MC-50E) Booklets with 25 Sheets per Booklet
- Approximate Sheet Dimensions: 4.9" x 2.9" (124 mm x 73 mm)



MC-50E



MC-5

These tissues can be used to clean optics or to protect an optic prior to placing it in a storage container. Tissues are free from contaminants and adhesives.

ITEM #	\$	£	€	RMB	DESCRIPTION
MC-5	\$ 9.80	£ 7.06	€ 8.53	¥ 78.11	Lens Tissue, 25 Sheets/Booklet, 5 Booklets
MC-50E	\$ 78.00	£ 56.16	€ 67,86	¥ 621.66	Lens Tissue, 25 Sheets/Booklet, 50 Booklets

Handi Pads

- Pure Cotton, Non-Woven, Lintless, Absorbent
- 100 Pads per Package
- Approximate Pad Dimensions: 4" x 8" (100 mm x 200 mm)



CP-100

These plush, non-woven cotton optical cleaning pads are free from contaminants and adhesives. Ideal for use with optic-cleaning solvents, they are Thorlabs' preferred cleaning pad.

ITEM #	\$	£	€	RMB
CP-100	\$ 11.40	£ 8.21	€ 9,92	¥ 90.86

Kimwipes

- Durable, Absorbent, Low-Lint Sheets
- 12 Boxes per Case, 280 Sheets per Box
- Approximate Sheet Dimensions: 4.5" x 8.4" (114 mm x 213 mm)



KW32

Kimwipes are extra-low-lint sheets commonly used for optic or fiber optic cleaning.

ITEM #	\$	£	€	RMB
KW32	\$ 44.60	£ 32.11	€ 38,80	¥ 355.46

Duster with Integrated Nozzle

Dusting gas is commonly used to clear particles off of an optic. This non-contact cleaning method is suitable for a variety of optics and consists of a trigger valve on a 10 oz can of Tetrafluoroethane (ozone safe).



CA3

Stocked in USA,
Germany, and Japan

Ground Shipment Required on Aerosol Cans

ITEM #	\$	£	€	RMB	DESCRIPTION
CA3	\$ 10.10	£ 7.27	€ 8,79	¥ 80.50	Duster with Plastic Nozzle

Premium Cotton-Tipped Applicators

- General Disposable Applicator
- Non-Sterile, Tightly Spun, Cotton Tip
- 6" Wood Stick with Cotton Tip
- Packaged in Quantities of 100



ITEM #	\$	£	€	RMB	DESCRIPTION
CTA10	\$ 4.20	£ 3.02	€ 3,65	¥ 33.47	100 Cotton-Tipped Applicators

Precision Optical Cleaner

This 3 oz ozone-safe can of optical-grade, fast-drying cleaning solution is designed to remove fingerprints, light oils, dust, fibers, water spots, and surface films from optical and plastic components. It is an ideal replacement for Isopropanol and Acetone.



FCS3

ITEM #	\$	£	€	RMB	DESCRIPTION
FCS3	\$ 16.40	£ 11.81	€ 14,27	¥ 130.71	Precision Optical Cleaner